

Imagista corporation



For support information
Tokyo Japan – October 14, 2025

Imagista corporation

Contact : Osamu Kobayashi

Tel/Fax 81-426-77-4745

kobayashi@imagista.co.jp

Revision information of X-ray image viewer XView_2025

Imagista corporation has announced revision of X-ray image viewer ver. 1.01 to new ver.2.01. In this revision, we further improved operability and added new features as described below. XView_2025 is image viewer software distributed to customers who purchase our products. It is designed to allow customers to quick and intuitive camera operation in settings camera parameter while viewing video image. Setting items are arranged on the same panel so results are quickly reflected in the image. This revision has the following points.

Add gamma correction and invert image movement checkbox

- Checking the checkbox will set the Gamma to default value 1.6.
- Checking the checkbox will invert the image vertically.



Organized functions (Save, Copy, 16BitSave, Exit)
They can be called from Menu Strip Item.

Add new functions Gamma and Upside Down